



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Multichannel Patient Monitoring System

Product Code :
SAMCOMEDICAL430001A002



Description :

Multichannel Patient Monitoring System We are engaged in manufacturing, supplying and exporting a wide range of Multichannel Patient Monitoring System. All our products are designed in compliance with international quality norms & standards, using premium raw material and latest technology. In addition to this, we check all the products on various parameters to assure flawlessness. Other information: 10.4" Color TFT display, with maximum 4 waveforms Parameters including ECG/RESP, SpO2, NIBP, TEMP and pulse rate Maximum 72-hour graphic and tabular trends of all parameters Suitable for adult, pediatric and neonatal patient Multichannel Patient Monitoring System, medical monitoring system exporters, medical monitoring system equipments manufacturers



Samco Medical

Website: www.samcomedical.com,