

**Product Name :**  
Syringe Infusion Pump

**Product Code :**  
SAMCOMEDICAL460001A002



#### Description :

Syringe Infusion Pump We are engaged in providing our clients with Syringe Infusion Pump that are available in different models and sizes at market leading prices. These products are widely appreciated by the customers for their corrosion & abrasion resistance, excellent performance, low maintenance and longer functional life. Furthermore, our products can be customized as per the specifications laid down by the clients. Product description: Self-explanatory with abbreviated operating instructions on LCD Easy-to-read liquid-crystal displays all essential data in large size numbers Reducing the chance of erroneous readings LCD provided with back light for viewing in dark environments Selectable three-step occlusion detection pressures Low, mid & high-occlusion setting "Keep-Vein-Open" function Infusion rate in KVO is 1% but not less than 0. 1 ml/hr This prevent the infusion line from blocking due to blood clots Preset Bolus Rate: A simple and safe Bolus function enabling the delivery of the shock does at any moment of the infusion with the preset rate as under, 300ml/hr for 10ml syringe, 500ml/hr for 20ml syringe, 700ml/hr for 30ml syringe, 1000ml/hr for 50/60ml syringe Full information about the course of infusion visual display of all the programmed parameters and infusion data such as syringe type and volume The infusion rate, the infused volume, time remaining till the end of infusion and occlusion pressure limit Ultimately 5 minutes prior the end of infusion Alarm announcing the necessity to prepare another dose of the drug Volume memory function In Memory Mode, this convenient feature retains the previously selected data of delivery rate and limit so they do have to be re-entered. This not only saves time and effort but also helps reduce errors External Switching mode power supply Complying with high standards, that assures a low and high voltage working environment, another important feature of the unit. A built in battery Guarantees safe operating of the pump in case of mains failure. When patient is being transported Syringe Infusion Pump, Infusion Pumps exporters, Infusion Pumps manufacturers

---

## Samco Medical

Website: [www.samcomedical.com](http://www.samcomedical.com),

