



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Vital Signs Monitor

Product Code :
SAMCOMEDICAL530001A001

Description :

Vital Signs Monitor 8" TFT, led backlight , 800x480 pixels. Visual and Audio alarms, adjustable. Patient management, name and age and id. Network with Central Station Software Audio volume adjustable from 1 ~ 8 level. Multi-language options: German, French, Italian, Turkish Built-in removable and rechargeable lithium battery UI with both Keypad and Touch Screen operation. Memory: Built-in memory or miniSD Memory card. Store more than 1000 pieces of archive Powerful data storage and management, Trend graph, Trend table, etc. Arrhythmia analysis and S-T segment analysis. ECG Lead mode: 3-lead or 5-lead HR range: 15 ~ 300 bpm Accuracy: $\pm 1/2$ bpm or $\pm 1/2$ %, which is greater NIBP Method: Oscillometry Mode: Manual/Auto/Continuous Measurement Range: 25 ~ 260mmHg SPO2 Measurement Range: 0 ~ 100% Accuracy: $\pm 1/2$ % (70%~100%) 0% ~ 69% unspecified Pulse Rate Measurement Range: 25 ~ 250 bpm Accuracy: $\pm 1/2$ bpm Dimension: 210 x 130 x 210mm Net weight: 2.3KG Vital Signs Monitor, Multipara Monitor Systems exporters, Multipara Monitor Systems manufacturers



Samco Medical

Website: www.samcomedical.com,