



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Digital Multipara Monitors

**Product Code :**  
SAMCOMEDICAL530001A003



**Description :**

Digital Multipara Monitors 8.4 = color TFT display Suitable for adult, pediatric and neonatal patient Basic parameters(ECG, SpO2, NIBP) in a durable case for bedside monitoring and transport monitoring Available and visual alarms with adjustable alarm ranges Net-workable data management and storage capacity ECG: Lead type:5-lead Input: RA; LA; RL; LL; V Sweep Speed: 2.5mm/s, 25mm/s, 50mm/s Accuracy:±1bpm or ±1% which is greater Protection: withstand 4000ac/50hz voltage in isolation against electro surgical and defibrillation Have S-T detection and Arrhythmia analysis facility Have Alarm, audible and visual alarm, alarm events re-callable facility Standard Configuration: ECG NIBP SpO2 Optional: RESP2-TEMP 2-IBP ET(02) Network Thermal Print Pacing Maker Digital Multipara Monitors, Multipara Monitor Systems exporters, Multipara Monitor Systems manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),