



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Fetal Monitor

**Product Code :**  
SAMCOMEDICAL520001A001

**Description :**

Fetal Monitors We are engaged in providing our clients with Fetal Monitors that is available in different models and sizes at market leading prices. All our products are widely appreciated by the customers for their corrosion & abrasion resistance, easy handling, low maintenance and longer functional life. Apart from this, our products can be customized as per the specifications laid down by the customers. Features: Wide 5.7 blue screen LCD Digital display and waveform oscilloscope of FHR and TOCO Continuous monitoring of Fetal Heart Rate and Uterine Contraction (OP) User friendly interface with menu function settings, direct and ease-of-use Preset usual functions shortcut keys for convenient operation Self-define multi-event marks Built-in high-resolution thermal array printer 7-chip impulse ultrasound transducer Highly sensitive Fetal movement auto detection Signal intensity indication ensures transducer detecting position. Auto power failure protection and reliable rechargeable High performance and self-adapting wide scope supply voltage design Full functional remote control Specifications: Uterine Contraction Pressure (TOCO) Measurement Measuring Method: External Measurement Measuring Range: 0~100 units Uterine Constriction pressure graphic display scope: 0~255 Uterine constriction pressure numeric display scope: 0~255 Auto Zero Function Yes. TOCO zero Points: 0, 5, 10, 15, 20, 25 & 30 Fetal Movement Record (Wireless): Detect Method: Manually Marked Wireless remote control switch Fetal Monitors, Fetal Monitors exporters, Fetal Monitors manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),