



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hand Held Single Channel ECG

Product Code :
SAMCOMEDICAL540001A001

Description :

Hand Held Single Channel ECG Input circuit: Floating Calibrating Voltage :1mV Input current: <50mA CMRR: >103dB A/D: 12bit Frequency Response: 0.05Hz--150Hz Filter: AC 50/60Hz EMG: 35-45Hz Skin Voltage Tolerance: +300mV Sensitivity Setting 5, 10, 20mm/mV, conversion deviation 25% Noisy Level: $\pm 15\mu V_{p-p}$ Safety standard: IEC class I, type CF Lead selection: Automatic and Manual Operation mode: Automatic and Manual Display: LCD 160*96 visual size: 34*60 Recording system: Thermal-array printing Power supply: AC220V $\pm 10\%$ DC 7.4 lithium battery Hand Held Single Channel ECG, ECG Machines exporters, ECG Machines manufacturers Size:190*90*40 Weight: About 800g



Samco Medical

Website: www.samcomedical.com,