



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
6 Channel ECG Machine

Product Code :
SAMCOMEDICAL540001A002



Description :

6 Channel ECG Machine
• Heart rate
• Heart rhythm
• For conduction abnormalities
• Prior heart attack
• Coronary artery disease
• Abnormally thickening of heart muscle
Features:
• 6 Channel ECG Machine are 7" TFT, touch screen, led backlight
• 6 Channel ECG Machine are Digital isolation technology and signal processing solution, digital filter
• 6 Channel ECG Machine are Patient management, name and age and id.
• 6 Channel ECG Machine are Detail analysis report.
• 6 Channel ECG Machine are Multi-language options: Germany, France, Italian, Turkey, etc
• 6 Channel ECG Machine are Built-in lithium rechargeable battery, wide thermal printing system
• 6 Channel ECG Machine are USB and LAN socket.
• 6 Channel ECG Machine are Memory: Built-in memory or mini SD Memory card. Store more than 1000 pieces of archive
• 6 Channel ECG Machine are Auto-measurement, auto-analysis and auto-interpretation
6 Channel ECG Machine, ECG Machines exporters, ECG Machines manufacturers



Samco Medical

Website: www.samcomedical.com,