



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
12 Channel ECG Machine

Product Code :
SAMCOMEDICAL540001A003



Description :

12 Channel ECG Machine Heart rate ½ Heart rhythm ½ For conduction abnormalities ½ Prior heart attack ½ Coronary artery disease ½ Abnormally thickening of heart muscle Feature: ½ Simultaneous acquisition of 12-lead ECG data, real-time recording and printing of 12-channel ECG waveform ½ Graphic display of 12-lead ECG waveform ½ Built-in analysis software of age which assures accurate analysis result ½ Auto-measurement, auto-interpretation, waveform playback and storage of ECG data ½ R5232 and USB interface cope with data-share or remote data management requirement ½ Option wireless function makes more convenient communication with PC ½ User friendly operation system contributed by built-in HELP function ½ Literal and graphic operation interface ½ Powerful filters to minimize interference ½ Heart rate measurement and pace-maker protection circuit ½ Multi-printing formats: manual/automatic, standard 12 channel, 3 channel plus 3 rhythm lead, 6 channel, 6 channel plus rhythm lead, 60s analysis of arrhythmia, R-R histogram, trend graph ½ AC, DC or built-in lithium battery power supply, alarm of battery weak and lead-off ½ Comfort and durability contributed by silicon gel keypad ½ Tremendous ECG data can be saved in built-in SD card ECG 330G 12 Channel ECG Machine, ECG Machines exporters, ECG Machines manufacturers



Samco Medical

Website: www.samcomedical.com,