



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Zoll PD 2000 Defibrillator

**Product Code :**  
SAMCOMEDICAL550001A001

#### Description :

Zoll PD 2000 Defibrillator Innovative features: Large Monitor Non-fade monitor ECG monitor Strip recorder AC power supply for stationary use Specifications: Type: VVI demand; asynchronous (fixed rate) when used without ECG leads Pulse Type: Rectilinear, constant current Pulse Duration: 40 milliseconds. Pulse Amplitude: Variable to 140 mA Pacing Rate: Variable from 30 to 180 PPM Output Protection: Fully defibrillator-protected and isolated Pacer On: Message display on monitor Defibrillator Waveform: Damped sinusoid Output Energy Delivered: Selectable at 2, 3, 5, 7, 10, 20, 30, 50, 100, 150, 200, 300, 360 joules Charge Time: Less than 10 seconds Delivered Energy Display: CRT monitor displays delivered energy Synchronized Mode: Synchronizes defibrillator pulse to patient's R-wave. SYNC message displayed on monitor. Marker on monitor and on recorder paper identifies R-wave Internal Paddles: Output limited to 50 joules when internal paddles are connected. Discharge from paddles Charge Controls: Control on apex paddle and front panel. Paddles: Standard paddles are anterior/anterior adult and pediatric. Adult paddles slide off to expose pediatric paddles Defibrillation Electrodes: Specifically designed ZOLL Multi- Function Electrodes used in anterior/anterior or anterior/posterior positions; packaged in pairs Integral Defibrillator Tester: Integral circuitry allows complete test of defibrillator charge and discharge without removing paddles from storage wells. Identical circuitry allows complete test of unit configured with multifunction cable Defib On: Message display on monitor Monitor and Display Patient Connection: Via 3-lead ECG cable, paddles, and electrodes. Selectable by front panel switch Input Protection: Fully defibrillator-protected. Special circuit prevents distortion of ECG by pacer pulse Bandwidth: .5 to 40 Hz (3dB) standard Display Format: Nonfade, moving trace with freeze capability Monitor On: Message display on monitor Screen Type: High-resolution CRT display Screen Size: 5 inches diagonally (108 mm x 85 mm viewing area) Zoll PD 2000 Defibrillator, Electrical Defibrillator exporters, Electrical Defibrillator manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),

