



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
HP Codemaster Defibrillator

Product Code :
SAMCOMEDICAL550001A002



Description :

HP Codemaster Defibrillator To fulfill the demands and needs of customers, we provide HP Codemaster Defibrillator that is available in professionally refurbished with a built-in recorder, ECG, pacing and AC power. Also, available with battery back-up, our products have several other functionalities. In addition to this, we provide customization as per the specifications laid down by the customers. Specifications: $\frac{1}{2}$ HP Codemaster Defibrillator are Dimensions: 7.9"H x 11.8"W x 15.6"L, 20cm x 30cm x 39.7cm $\frac{1}{2}$ HP Codemaster Defibrillator are Weight: 24 lbs. / 10.9kg Battery: $\frac{1}{2}$ Type: Rechargeable, sealed lead-acid battery. 4 Ah, 12 V nominal. $\frac{1}{2}$ Charge Time: 2 hours to 90% of full capacity. 18 hours to 100% capacity. Repeated charging to less than 100% will reduce the useful life of the battery. $\frac{1}{2}$ Capacity: 2.5 hours monitoring or 50 full-energy discharges or 1 hour monitoring and recording. $\frac{1}{2}$ AC input: 100 to 230 VAC +/- 15%, 50/60 Hz. $\frac{1}{2}$ Battery Indicators: Illuminated LED indicates battery is charging. Low Battery message appears on the monitor Monitor: $\frac{1}{2}$ Input: ECG activity may be viewed through the paddles or patient cable $\frac{1}{2}$ Lead fault: Leads off message and dashed baseline appear on the monitor $\frac{1}{2}$ Common mode rejection: Greater than 100dB measured as per AAMI standards for cardiac monitors HP Codemaster Defibrillator, Electrical Defibrillator exporters, Electrical Defibrillator manufacturers



Samco Medical

Website: www.samcomedical.com,