

Product Name :
Versamed Ivent 201 Ventilator

Product Code :
SAMCOMEDICAL560001A008



Description :

Versamed Ivent 201 Ventilator Electrically Powered, Computer Controlled Ventilator Closed Loop Flow & Pressure Controller Ventilation Modes Assist Control (A/C) Volume Controlled Pressure Controlled Synchronized Intermittent Mandatory Ventilation (SIMV): Volume Controlled Pressure Controlled Continuous Positive Airway Pressure (CPAP) Pressure Support Ventilation (PSV) Special Features: Adaptive Bi-Level (NIV) Adaptive Flow^{1/2} Adjustable Rise Time Adaptive Time^{1/2} Easy Exhale^{1/2} Preset Parameters by Patient Weight Sigh Breath Apnea Backup Ventilation 100% O₂ Suction Ventilation Performance and Controlled Parameters: Respiratory Rate : 1 to 80 BPM Tidal Volume : 50 to 2,000 ml Inspiratory Pressure Limit : 5 to 80 cm H₂O Inspiratory Time Adaptive Time^{1/2} : , or 0.2 to 3 seconds Peak Flow Adaptive Flow^{1/2} : or 1 to 120L/min Spontaneous up to 180 L/min Oxygen Mix (FIO₂) : 21% to 100% PEEP : 0 to 40 cm H₂O Trigger Sensitivity : 1 to 20 L/min Flow Sensitivity + off ^{1/2}0.5 to ^{1/2}20 cm H₂O Pressure Sensitivity + off PSV : 0 to 60 cm H₂O Positive Pressure Relief Valve : 80 cm H₂O Synchronized Nebulizer: 10-240 Minutes Monitoring and Displayed Parameters: Airway Pressure (analog bar graph & numerical) Total Breath Rate I:E Ratio Exhaled Tidal Volume Exhaled Minute Volume Peak Flow Inspiratory Time Electrical Power Source (external / internal) Battery Level Pressure, Flow and Volume Waveforms Software Package: Real Time Pressure and Flow Waveforms Waveform History Browse Trending of Monitored Data (72 Hours) Respiratory Mechanics (C, R, MAP, RR/Vt) Pressure, Flow and Volume Loops Lung Mechanics (C, R, MAP, RR/Vt) Versamed Ivent 201 Ventilator, Medical Ventilator exporters, Medical Ventilator manufacturers

Website: www.samcomedical.com,

