



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
LP10 Ventilator

**Product Code :**  
SAMCOMEDICAL560001A009



**Description :**

LP10 Ventilator  
• Pressure control knob for pressure-limited, time-cycled ventilation  
• Decreased internal resistance reduces breathing effort in SIMV mode  
• Continuous, low-maintenance operation with minimal supervision results in low cost of ownership  
• Features and portability provide flexibility to meet a broad range of ventilation and patient requirements  
• Internal and external battery operation provides convenience and security  
Indications For Usage:  
• A portable ventilator ideal for in hospital transport, sub acute facilities and long term, in-home care



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),