



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Drager Savina

**Product Code :**  
SAMCOMEDICAL560001A011



**Description :**

Drager Savina Ventilation modes - IPPV (CMV), IPPV Assist (CMV Assist)  $\dot{V}_T$  SIMV, SIMVASB (SIMV/PS)  $\dot{V}_T$  CPAP, CPAPASB (CPAP/PS)  $\dot{V}_T$  BIPAP1), (PCV+) (optional), BIPAP1)ASB (PCV+/PS) (optional)  $\dot{V}_T$  Enhancements - NIV  $\dot{V}_T$  Non Invasive Ventilation with optimized alarm system and automatic leakage compensation (optional)  $\dot{V}_T$  Auto Flow  $\dot{V}_T$  - Automatic adaptation of the inspiratory flow in volume orientated ventilation modes (optional)  $\dot{V}_T$  LPO - Low Pressure Oxygen. Independent oxygen supply, e.g. with an O2 concentrator (optional)  $\dot{V}_T$  Graphic screen - Advanced ventilation monitoring (optional)  $\dot{V}_T$  Nurse call - Connection for transmitting alarm signals to a central Drager Savina, Medical Ventilator exporters, Medical Ventilator manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),