



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Pulse Oximeter System

**Product Code :**  
SAMCOMEDICAL570001A001

**Description :**

Pulse Oximeter System Family Hospital Oxygen bar Social medical organization Physical care in sports  
Features: Large screen: SpO2,PR are shown on LED screen, Pleth waveform is displayed on 3.5" LCD screen in real time Light in weight, only 1kg with battery 72-hour data storage, 100-item patient ID 72-hour SpO2 and PR spot graphs 3-level visual & audio alarm Intellectual control of power consumption, internal rechargeable battery lasts 24 hours or AC power Battery-low indication Built-in thermal recorder(optional) Suitable for adult, pediatric and neonatal Widely used in small clinic, all departments in hospital, ambulance Compatible with CHOICE central monitoring system USB data upload (optional)



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),