



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Syringe Pump Dual Channel

Product Code :
SAMCOMEDICAL590001A001

Description :

Syringe Pump Dual Channel International standard syringe pump Compatible: Automatic syringe calibration, suitable for all brands Convenient: $\frac{1}{2}$ Openable slider hook $\frac{1}{2}$, convenient to clamp syringe plunger Powerful: Battery capable to support more than 6 hours for one channel continuous working Energy saving: Idle channel can be set to sleeping mode Function: Dual-channel with different flow rate Channel sleeping function 3-level occlusion alarm limits setting Limits setting Automatically reduce impact dose release after occlusion Delivered volume inquiry KVO flow rate Fast rate control $\frac{1}{2}$ Purge $\frac{1}{2}$ button activate safe protection Flow rate setting button lock Abundant optic and sound alarm AC and internal power indication Standard RS232 interface 500 history records Syringe Pump Dual Channel, Syringe Pump exporters, Syringe Pump manufacturers



Samco Medical

Website: www.samcomedical.com,