



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Radiant Warmer with Detachable Trolley

Product Code :
SAMCOMEDICAL600001A001

Description :

Radiant Warmer with Detachable Trolley We are counted among the reputed organizations, engaged in offering our clients with Radiant Warmer with Detachable Trolley. All our products are available in different sizes, shapes and designs at market leading prices to fulfill the demands and needs of customers. Apart from this, we stringently check all the products on various parameters to ensure flawlessness. Features: $\frac{1}{2}$ Easy usage $\frac{1}{2}$ Low maintenance $\frac{1}{2}$ Lightweight Mechanical Specifications $\frac{1}{2}$ Easy access of Infant from all sides with thick Acrylic collapsible side panels $\frac{1}{2}$ Acrylic baby tray with Infant's Head up / down facility & foam mattress. $\frac{1}{2}$ Provision for positioning X - ray cassette. Technical Specifications $\frac{1}{2}$ Power Source: 60 watts bulb to examine babies $\frac{1}{2}$ Light Source: 600 watts $\frac{1}{2}$ Heater Capacity: 800 watts at 100% Heater power $\frac{1}{2}$ Power Consumption: 230V AC, $\frac{1}{2}$ 10%, 50Hz



Samco Medical

Website: www.samcomedical.com,