

Product Name :
LED Phototherapy Stand

Product Code :
SAMCOMEDICAL600001A002



Description :

LED Phototherapy Stand
 LED Phototherapy Stand
 Easy to operate
 Low power consumption (12 W) compared to any conventional phototherapy system (120 W)
 LED's emitting blue light in 450-470 nm. Spectrum matching peak absorption wavelength (458 nm.) at which bilirubin is broken down
 Do not emit light in infrared radiation range reduction the risk of fluid loss
 Long life > 20,000 hrs
 Negligible Maintenance cost
 Intensity: > 25 W / cm² / nm
 Heavy duty mobile stand with hight adjustment facility & easy swivelling box when stand is to be used in conjunction with radiant warmer or infant trolley
 No need to Cover eyes of the baby
 Mounted on casters for easy mobility
 M. S. Fabrication with powder coated paint
 Time Totalizer (Optional)
 Operation Voltage: 100-250V AC, 50Hz
 Contact Phototherapy:
 Compact, Portable & Easy to Carry
 Detachable Baby Tray for Cleaning
 No need to Cover eyes of the baby
 Low Power consumption (12 W)
 Led's emitting blue light in 450-470 nm. Spectrum matching peak absorption wavelength (458 nm.) at which bilirubin is broken down
 Do not emit light in infrared radiation range reduction the risk of fluid loss
 Long life > 20,000 hrs
 Irradiance: > 25 W / cm² / nm
 Can be used separately in mother bed / in trolley
 M. S. Fabrication with powder coated paint
 Time Totalizer (Optional)
 Operation Voltage : 100-250V AC, 50Hz
 LED Phototherapy Stand, Neonatal Products exporters, Neonatal Products manufacturers



Samco Medical

Website: www.samcomedical.com,

