



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
X Ray Machine

Product Code :
SAMCOMEDICAL610001A003



Description :

X Ray Machine Single phase Full wave rectified (silicon rectifiers) H medium powered X-Ray unit independent KV mA and time techniques for better radio graphic results Booster Transformer for constant mA output Instantaneous electronic tube overload protection for longer tube life Independent KV, mA and time techniques for better radio graphic results Voltage compensation with coarse and fine voltage compensation switch 6 LED indicators for Line, X-Ray, 0/L, bucky, standby and fuse Precision solid state electronic timer Automatic calculation of mAS with digital display Digital display of KVP



Samco Medical

Website: www.samcomedical.com,