



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
100MA Mobile X Ray Machines

Product Code :
SAMCOMEDICAL610001A004



Description :

100MA Mobile X Ray Machines Counter-balanced. tube arm for increased reach better positioning Integrated cassette box for easier operation and better cassette handling Elegant digital display and extremely easy to use controls MA meter & voltmeter for line and exposure monitoring Automatic calculation of mass with digital display Digital display of KVP 6 Led indicator for line, 0/L, X Ray, Stand by, Bucky & fuse Precision solid state electronic time.



Samco Medical

Website: www.samcomedical.com,