



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
300ma X Ray Machines

Product Code :
SAMCOMEDICAL610001A005



Description :

300ma X Ray Machines Radio graphic output: 40kw as per IS7620 (Parts) Tube filament current stabilization through constant voltage transformer (CVT) Full wave rectification Digital display of KV, mA and MAS Overload protection (Blockage of Exposure) Shortest exposure time 10 MAS in 500mA Supports serialography up to 6fps kv, mA and MAS values set with step switches fluoroscopy time totalizer 5 LEO indicators for Line, 0/L exposure, ready and bucky in voltage stabilization (compensation) Through course and fine voltage compensation switches voltmeter and mA meter for viewing input line voltage and mA during exposure 300ma X Ray Machines, X-Ray Machine & Instrument exporters, X-Ray Machine & Instrument manufacturers



Samco Medical

Website: www.samcomedical.com,