



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Corner Connected X-Ray Cassettes

Product Code :
SAMCOMEDICAL610001A006



Description :

Corner Connected X-Ray Cassettes Medical X-Ray Cassettes for Radiography Available in different sizes. Unit inches & Cms The Cassettes are designed considering the latest technology which results with improved image quality The body frame is made up of Aluminium for light weight connected on-four sides with non metallic construction to give strong cover soft push buttons locking With stainless steel springs & high grade foam laminated placed on back board for quality results of film contrast which can be easily replaceable The lead sheet on the back board to eliminate phenomena of X-ray film



Samco Medical

Website: www.samcomedical.com,