



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
X-Ray Film Viewing Boxes

Product Code :
SAMCOMEDICAL610001A007



Description :

X-Ray Film Viewing Boxes The view Boxes designed to reduce glare and excessive heat provide with uniform illumination Other feature include ventilation for cooling, a drip tray for wet films. Available in single, double and four film viewings boxes Sizes: LDH 15 X 5 X 25 inch (Single film Viewer) 24 11/2 5 X 19 1/2 inch (Double Film Viewer) 48 1/2 X 5 X 19 1/2 inch (Four Film Viewer) X-Ray Film Viewing Boxes, X-Ray Machine & Instrument exporters, X-Ray Machine & Instrument manufacturers



Samco Medical

Website: www.samcomedical.com,