



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
60mA 100mA X-Ray Equipments

**Product Code :**  
SAMCOMEDICAL610001A008



**Description :**

60mA 100mA X-Ray Equipments Eco Mobile-Fully counter balanced very compact-space saver Articulated Tube Head Arms to position the Tube Head for all possible radiographic positions. Fully counter balanced Tube Head for ease of handling with 360° rotation for column Floor to Ceiling column with straight cross arm - Ideal for radiographic work Floor to ceiling column with  $\frac{1}{2}$ U $\frac{1}{2}$  arm provides screening & with SFD provides spot filming One position diagnostic table with / without Bucky Five position diagnostic table having 00, 300, 600, 900 & - 150, trendelenberg position Salient Features : Fullwave rectification which delivers almost double the output than the self rectified X - Ray Systems. Full wave rectification with high voltage silicon rectifier for low losses and consistent radiographic results Smoothing circuit provides maximum mA output accuracy at all time settings Full field light beam Collimator - for proper positioning and preventing accidental cut off the field Continuous fluoroscopic timer by pressing F\_RST push button. (100mA) Low output power station 4Kw, for versatile, in absence of power-Line. (100mA) Double mA range meter with lamp. (100mA). 60mA 100mA X-Ray Equipments, X-Ray Machine & Instrument exporters, X-Ray Machine & Instrument manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),