



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Single Beam Spectrophotometer

Product Code :
SAMCOMEDICAL650001A006



Description :

Single Beam Spectrophotometer Optical System Single Beam, Grating 1200 lines/mm Wavelength Range 200-1000nm Band Width 1 nm Wavelength Accuracy ± 0.5 nm Wavelength Repeatability ± 0.1 nm Wavelength Setting Automatic Photometric Accuracy ± 0.1 %T Photometric Repeatability ± 0.3 %T Photometric Display Range -0.3-3.0A, 0-200%T, 0-9999C Stability ± 0.002 A/h @ 500nm Stray Light ≤ 0.3 %T @ 220nm, 360nm Data Output Port USB Printer Port Parallel Port Display 128*64 Dots LCD Lamps Deuterium Lamp & Tungsten Halogen Lamp Detector Silicon Photodiode Power Requirement AC220V/50 Hz Single Beam Spectrophotometer, Double Beam UV-VIS Spectrophotometer exporters, Double Beam UV-VIS Spectrophotometer manufacturers



Samco Medical

Website: www.samcomedical.com,