



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
UV Spectrophotometer

Product Code :
SAMCOMEDICAL660001A003



Description :

UV Spectrophotometer Single Beam Microprocessor UV/VIS Spectrophotometer Range 195 to 1000 nm with facility for automatic concentration, % Transmission, Absorbance & K Factor measurement. It has sample holder with 4 position selector control for 10 mm path length cuvettes. Facility has been provided for interface with any printer (RST 232C). Set of 2 quartz and 4 glass cuvettes is provided. SOFTWARE INCLUDED (Imported) UV Spectrophotometer, Single Beam UV-Vis Spectrophotometer exporters, Single Beam UV-Vis Spectrophotometer manufacturers



Samco Medical

Website: www.samcomedical.com,