



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Atomic Absorption Spectrophotometers

Product Code :
SAMCOMEDICAL690001A001

Description :

Atomic Absorption Spectrophotometer ½ Automatically controlled changeover of the Integrated flame and graphite furnace automizer featuring easy operation and time saving eliminated human labor ½ A flame emission burner head can be installed to perform flame emission analysis to Alkali metals as K, Na etc. Accurate Fully Automated Control System: ½ Automatic 6 lamp turret, automatic adjustment of lamp current and optimization of light beam position ½ Automatic wavelength scanning and peak picking ½ Automatic special spectral bandwidth changing ½ Automatic changeover between flame and graphite furnace operation, automatic optimization of position parameters automatic ignition and automatic gas flow setting Atomic Absorption Spectrophotometer, Atomic Absorption Spectrophotometer exporters, Atomic Absorption Spectrophotometer manufacturers



Samco Medical

Website: www.samcomedical.com,