



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Digital pH Meters

Product Code :
SAMCOMEDICAL700001A001

Description :

Digital pH Meters Display 3 1/2 digit LED 3 1/2 digit LED 3 1/2 digit LCD Range pH 0 to 14.00 0 to 14.00 0 to 14.00 mV 0 to + 1999 0 to + 1999 - Resolution pH 0.01 0.01 0.01 mV 1 1 - Accuracy pH +0.01 +0.01 +0.01 mV +1 +1 - Temperature Auto 0 to 100o - - Compensation Manual 0 to 100o Manual 0 to 100o Manual 0 to 100o Slope Correction 80 to 120% 80 to 120% 80 to 120% Power 230V + 10% 230V + 10% 1 x 9V Battery Dimensions 76 x 275 x 175 76 x 275 x 175 85 x 28 x 165 Accessories: 1/2 Combination pH Electrode 1/2 Buffer Tablets 4pH & 7pH 1/2 Operation Manual 1/2 Electrode Stand 1/2 Dustcover Digital pH Meters, PH Meters exporters, PH Meters manufacturers



Samco Medical

Website: www.samcomedical.com,