



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Digital Hemoglobin Meter

**Product Code :**  
SAMCOMEDICAL730001A002



**Description :**

Digital Hemoglobin Meter Instrument Digital Hemoglobin Meter Display Digital 3 Digit, Bright LED Display Mode % T, O.D. and 3 concentration Modes, C1,C2,C3 Zero Adjustment Automatic (Without Knob rotation) Accuracy + 0.02 O.D. Stability + 0.02 O.D. Per Hour Light Source 6.8 Volt, 0.3 Amp Tungsten Lamp Amplifier Micro controller Based Memory In built memory to retain last K factor for concentration even after power loss. Sample Volume 1 ml. minimum Power 230 V + 10 % AC, 50 Hz. Accessories Set of Test Tube- 5 No. Dust Cover 1 No. Operation Manual Salient Features: Fully Automatic Digital Read Out Highly Accurate Long Life Digital Hemoglobin Meter, Hemoglobin Meter exporters, Hemoglobin Meter manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),