



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Flame Photometer

Product Code :
SAMCOMEDICAL750001A001

Description :

Flame Photometer Digital Flame Photometer With fitted (Na & K) 3½ digit LED display (ppm) along with compressor and other accessories Ca & Li filter can be provided on request Digital Clinical Flame Photometer Dual channel single aspiration with filter (Na & k) for medical use. Final result in meq/L 2½ digit dual LED display for simultaneous results along with compressor & accessories Microprocessor flame photometer With Sodium (Na) and Potassium (K) filters, upto 5 points calibration, direct readout in pp, & meq/l, auto filter selection, auto ignition, 24 x 4 line LCD display, Single Aspiration, RS232 interface, printer attachments facility, with compressor and other accessories for lab & medical use. Ca & Li Filters and optional at extra cost Flame Photometer, Digital Flame Photometers exporters, Digital Flame Photometers manufacturers



Samco Medical

Website: www.samcomedical.com,