



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Digital Friability Tests Apparatus

**Product Code :**  
SAMCOMEDICAL770001A001

**Description :**

Digital Friability Test Apparatus Digital Friability Test Apparatus (Double Drum) Measures Obsession Strength of Tablets with programmable timer. Rotation Speed 25 rpm +/- 1, auto & manual mode of operation with 2 acrylic drums, Results are displayed on 3 digit LED display and it has 3 digit prestable control Digital Friability Test Apparatus (Single Drum) Measures Obsession Strength of Tablets with programmable timer. Rotation Speed 25 rpm +/- 1, auto & manual mode of operation with 1 acrylic drums, Results are displayed on 3 digit LED display and it has 3 digit prestable control Digital Friability Test Apparatus, Digital Friability Test Apparatus exporters, Digital Friability Test Apparatus manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),