



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Digital Conductivity Meters

**Product Code :**  
SAMCOMEDICAL780001A001

**Description :**

Digital Conductivity Meters Display 3 1/2 Digit LED 3 1/2 Digit LCD Range 5 Ranges 0 - 1000 Mohs/cm 2 Ranges 2 & 20 Mohs/cm Temperature Compensation Manual 0-50oC - Cell Constraint Adjustment Adjustable Pre Adjustable Conductivity Cell Platinum DIP Platinum DIP Measuring Frequency 1000 Hz 1000 Hz Accuracy +1% FSI +1% FSI Power 230v + 10% 1 x 9V Battery Dimensions 76 x 275 x 175 mm 85 x 28 x 165 mm Weight 2 kg. (Approx.) 0.8 kg. (Approx.) Accessories Conductivity Cell, Operation Manual, Dust Cover Digital Conductivity Meters, Digital Conductivity TDS Meters exporters, Digital Conductivity TDS Meters manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),