



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Microprocessor Conductivity TDS Meter

**Product Code :**  
SAMCOMEDICAL780001A002



**Description :**

Microprocessor Conductivity TDS Meter  $\frac{1}{2}$  0 to 1000 milli mho in 5 ranges Range Selection Automatic:  $\frac{1}{2}$  0 to 20 micro mho  $\frac{1}{2}$  20 to 200 micro mho  $\frac{1}{2}$  200 to 2 milli mho  $\frac{1}{2}$  2 to 20 milli mho  $\frac{1}{2}$  20 to 200 milli mho  $\frac{1}{2}$  200 to 1000 milli mho  $\frac{1}{2}$  Resolution 0.1 micro mho  $\frac{1}{2}$  Cell Constant Adjustable  $\frac{1}{2}$  Temperature co-eff. Adjustable between 0 to 4%  $\frac{1}{2}$  Temperature Comp. Automatic for entire range  $\frac{1}{2}$  Accuracy +0.5% of range + 1 digit TDS:  $\frac{1}{2}$  Range : 0 to 200 ppt in 5 ranges  $\frac{1}{2}$  Range Selection : Automatic  $\frac{1}{2}$  Resolution : 0.1 ppm  $\frac{1}{2}$  TDS Factor : 0.5 and adjustable manually Salinity:  $\frac{1}{2}$  Range : 0 to 50 ppt  $\frac{1}{2}$  Resolution : 0.1 ppt  $\frac{1}{2}$  Accuracy : +0.5% of range + 1 digit  $\frac{1}{2}$  Temperature Comp. : Automatic  $\frac{1}{2}$  Salinity Factor : .75 and Adjustable Manually Temperature:  $\frac{1}{2}$  Range -200 TO +200  $\frac{1}{2}$ C  $\frac{1}{2}$  Resolution 0.1  $\frac{1}{2}$ C  $\frac{1}{2}$  Accuracy + 0.2 %  $\frac{1}{2}$  Input Type Pt 100  $\frac{1}{2}$  Calibration Manual  $\frac{1}{2}$  Display : 16 Character \* 2 Lines Alphanumeric Lcd.  $\frac{1}{2}$  Display Format : Temp & Cond Or Temp & Tds, Temp & Sal Simultaneously  $\frac{1}{2}$  Setting : By Soft Touch Keys  $\frac{1}{2}$  Storage : 1000 Samples With Date & Time  $\frac{1}{2}$  Interface : Rs-232, Centronics & Computer  $\frac{1}{2}$  Real Time Display : 24 Hour Mode With Date  $\frac{1}{2}$  Error Indication : T- Error, Cal Error  $\frac{1}{2}$  Power Supply : 230 V  $\frac{1}{2}$  10% Ac 50 Hz  $\frac{1}{2}$  Digital Output : A Serial Rs-232 B Parallel Centronics  $\frac{1}{2}$  Print Format : A All Data B Date Wise Microprocessor Conductivity TDS Meter, Digital Conductivity TDS Meters exporters, Digital Conductivity TDS Meters manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),

