



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dissolved Oxygen Meters

Product Code :
SAMCOMEDICAL790001A001

Description :

Dissolved Oxygen Meters Range D.O. 0 to 20.0 ppm Resolution D.O. 0.1 ppm Accuracy D.O. 0.1 ppm + 1 count Temp. Compensation 0o C to 50o C Manual Compensation D.O. 3.00% D.O. Sensor type D.O. Amperometric Gold/Silver membrane General: Display 3 1/2 digit LED display Power 230V A.C. 50 Hz Weight 2 Kg Standard Accessories: DO Electrode 1 Nos. Instruction Manual 1 Nos. D.O. Membranes 10 Nos. Dissolved Oxygen Meters, Digital D.O. Meters exporters, Digital D.O. Meters manufacturers



Samco Medical

Website: www.samcomedical.com,