



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Microcontroller pH Meters

Product Code :
SAMCOMEDICAL820001A001

Description :

Microcontroller pH Meter Range 0 TO 14 pH Resolution 0.01 pH Accuracy + 0.01 pH + 1 digit Slope Adjustment 80 % to 120 % by default 100 %. Temp. Compensation 0 to 100 °C Auto with temp. probe/Manual with Keys. Auto Buffer Recognition 4.00 pH, 7.00pH, 9.20 pH (AUTOMATIC) Calibration AUTO & MANUAL, 3 Point/2 Point (Programmable) 2- POINT, 3- POINT mV (ORP) Range - 1999 mV to + 1999 mV Resolution 1 mV Accuracy + 01% + 1 digit Temperature Range -50.0 TO +200 °C Resolution 0.1 °C Accuracy + 0.1 % of range + 1 digit Input Type Pt 100 probe Data Entry By Six soft touch keys with audible sound Storage 1000 Samples Real time display 24 Hour mode with Date Interface RS 232 & Parallel Print Format All Data & Date Wise (Programmable) Digital Output Serial RS 232 & Printer Centronics (Programmable) Microcontroller pH Meter, Microcontroller pH Meter exporters, Microcontroller pH Meter manufacturers



Samco Medical

Website: www.samcomedical.com,