



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Air-Compressing Nebulizers

Product Code :
SAMCOMEDICAL830001A001

Description :

Air-Compressing Nebulizers This equipment sprays liquid medication to fog panel by compressing air, and forms tiny particles, which flow into the throat through imbibing tube. It adopts low noise, oil-free high-efficiency valve with no need for daily lubrication. Ideal product for family and medical units' use & suitable for all kinds of patients. $\frac{1}{2}$ Power supply: AC220 $\frac{1}{2}$ 20V, 50 $\frac{1}{2}$ 1HZ $\frac{1}{2}$ Max negative pressure compressing pump: ? 0.15MPa $\frac{1}{2}$ Free air flux compressing pump: ?10L/min $\frac{1}{2}$ Max nebulizing rate: ?0.10mL/min $\frac{1}{2}$ Noise: ?65dBA



Samco Medical

Website: www.samcomedical.com,