



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
32 Channel EEG

**Product Code :**  
SAMCOMEDICAL850001A002



**Description :**

32 Channel EEG • 32 Channel simultaneous acquisition of raw data • User programmable acquisition and photic sequences protocols • Online true AC impedance check • Automatic paging facility for quick review and analysis • Brain mapping color coding as per international standard • EEG synchronous Video Recording (optional) • Available in portable model also 32 Channel EEG, Neurology Devices exporters, Neurology Devices manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),