



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Fixed Machine with Rotating Anode Tube

**Product Code :**

SAMCOMEDICAL850001A003

**Description :**

Fixed Machine with Rotating Anode Tube 1/2 100mA, 100 KVP, X-Ray Machine for Radiography only/ Rad + Fluoro applications 1/2 8 KW X-Ray Generator 1/2 Electronic Timer 1/2 Electronic Overload Protection Device 1/2 Collimator 1/2 Floor to ceiling Stand 1/2 Tubes: Rotating Anode X-Ray Tube : BEL or Toshiba Tables: 1/2 Horizontal Bucky Table (HRT) 1/2 Horizontal Table with Floating Top (HFT) 1/2 5 position Manual Tilt Table (MPT) Fixed Machine with Rotating Anode Tube, Neurology Devices exporters, Neurology Devices manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),