



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Mobile X-Ray Machine with Anatomical Programs

Product Code :

SAMCOMEDICAL860001A001

Description :

Mobile X-Ray Machine with Anatomical Programs $\frac{1}{2}$ 100mA, 100 KVP, X-Ray Generator $\frac{1}{2}$ 1st in India of its Kind. Microcontroller Based Design $\frac{1}{2}$ Soft Touch control Panel $\frac{1}{2}$ Graphical LCD Display $\frac{1}{2}$ Anatomical programme (120) $\frac{1}{2}$ Electronic Overload Protection, Error Indication $\frac{1}{2}$ Stationary Anode X-Ray Tube : BEL or Imported $\frac{1}{2}$ Counter Balanced Mobile Stand $\frac{1}{2}$ Suitable for Bed Side X-rays Mobile X-Ray Machine with Anatomical Programs, Radiology Devices exporters, Radiology Devices manufacturers



Samco Medical

Website: www.samcomedical.com,