



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
100 mA, 100 KVP Fixed X-Ray Machine

**Product Code :**  
SAMCOMEDICAL860001A003



**Description :**

100 mA, 100 KVP Fixed X-Ray Machine ÷ 100mA, 100 KVP, Fixed X-Ray Machine for Radiography & Fluoroscopy ÷ 8 KW X-Ray Generator ÷ Electronic Timer ÷ Electronic Overload Protection Device ÷ Stationary Anode X-Ray Tube : BEL or Imported ÷ Collimator ÷ Tubes stand:- Floor to ceiling Tables: ÷ orizontal Bucky Table (HRT) ÷ Horizontal Table with Floating Top (HFT) ÷ 5 position Manual Tilt Table (MPT) 100 mA, 100 KVP Fixed X-Ray Machine, Radiology Devices exporters, Radiology Devices manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),