



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
300 MA, 125 KVP X-Ray Machine

**Product Code :**  
SAMCOMEDICAL860001A004



**Description :**

300 MA, 125 KVP X-Ray Machine  
Applications  
24 KW X-Ray Generator  
Electronic Timer  
Electronic Overload Protection Device  
Floor to ceiling Stand  
Tubes: Rotating Anode X-Ray Tube : BEL or Toshiba  
Tables: Horizontal Bucky Table (HRT)  
Horizontal Table with Floating Top (HFT)  
5 position Manual Tilt Table (MPT)  
Motor Driven table (MDT)  
300 MA, 125 KVP X-Ray Machine, Radiology Devices exporters, Radiology Devices manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),