



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
500 MA, 125 KVP X-Ray Machine

Product Code :
SAMCOMEDICAL860001A005



Description :

500 MA, 125 KVP X-Ray Machine
• 500mA, 125 KVP X-Ray Machine for Radiography Only or Rad + Fluoro Applications
• 40 KW X-Ray Generator
• Electronic Timer
• Electronic Overload Protection Device
• Floor to ceiling Stand
• Tubes: Rotating Anode X-Ray Tube : BEL or Toshiba
• Tables: Horizontal Bucky Table (HRT)
• Horizontal Table with Floating Top (HFT)
• 5 position Manual Tilt Table (MPT)
• Motor Driven table (MDT)
500 MA, 125 KVP X-Ray Machine, Radiology Devices exporters, Radiology Devices manufacturers



Samco Medical

Website: www.samcomedical.com,