



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Post Mortem knife

Product Code :
SAMCOMEDICAL940001A001

Description :

Post Mortem Knives Blade numbers PM 40 & PM 40B are available in heavy duty Carbon Steel. Supplied non Sterile, individually wrapped in a V.C.I. paper. 10 pieces in a box. DISPOSABLE POST MORTEM KNIVES Disposable Post-Mortem handles are available in blade nos. PM 40 & PM 40B. BLADE NO. 60 & 70 Available in both carbon & stainless steel supplies sterile & non-sterile in a box of 100 pcs. POST MORTEM KNIFE HANDLES Post-Mortem Handle take blade nos. PM 40 & PM 40B, Handle No. 8 fit blade No. 60 & No. 70.



Samco Medical

Website: www.samcomedical.com,