



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Blood And Surgical Fluid Absorbent Floor Pads

Product Code :
SAMCOMEDICAL950001A001

Description :

Blood And Surgical Fluid Absorbent Floor Pads We offer Blood and Surgical Fluid Absorbent Floor Pads that can reduce clean-up time and protect health care workers against exposure to infectious waste. The olefin non-skid backing of these floor pads can easily lower the instances of slippages and prevents the pad from curling at the edges. Moreover, this olefin non-skid backing keeps floor under pad dry and stain free and prevents bunching under stool wheels. The features of absorbent non-woven material of these pads are as follows : Optimum wicking and retention of fluids. Low linting Light weight reduces disposal cost Safe for incineration Clients can avail these products from us in four convenient sizes : 36"x30" ½ 25 per box 36"x33" ½ 25 per box 36"x40" ½ 25 per box 30"x72" ½ 25 per box Blood And Surgical Fluid Absorbent Floor Pads, Consumable Operation Theatre Items exporters, Consumable Operation Theatre Items manufacturers



Samco Medical

Website: www.samcomedical.com,