



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Radiant Warmer, Open Care System

Product Code :
SAMCOMEDICAL130001A003

Description :

Infant Radiant Warmer, Open Care System

Technical Specification :

Infant Radiant Warmer (Infant Resuscitation Trolley), Open Care System on wheels. Consists of Ceramic Infra Red Heater 650 Watts. Metallic overhead heater box. The construction of the stainless steel reflector and guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards. Halogen Examination Lamp (12 Volts, 50 Watts). To provide accurate assessment of Infant colour and sufficient illumination during procedures. Essential Attachments: S.S. IV Stand. Infant Care Trolley Standard. Clear acrylic collapsible side panels to assess an infant from all sides. Stainless steel baby tray with transparent slot to incorporate the facility of under surface photo-therapy. Acrylic tray and cushion mattress with head up/down facility. Acrylic X-ray tray with X-ray guide for positioning the cassette. Two Lower shelves. Mounted on wheel 2inch, The design aspects also allow easy of dismantling baby care area during cleaning or fumigation. Feature: Control Modes: Skin / Air / Manual Range: Skin control: 34 to 38 C Air control: Room temp. to 39 C Resolution: 0.1 C Accuracy: ± 0.2 C Display: 1in bright skin & air temp.(7 segment LED) Status: On LCD Control Mode & Set temp Timer: Timer / APGAR Timer / Baby Attend Manual mode: Time Programmable Heater output settable Over temp. cut-off Alarm: Auto audio alarm mute, Mute time programmable Safety alarms: ON LCD Heater Failure Power Control Failure Probe & Power failure alarm Infant Radiant Warmer, Open Care System, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,