

Product Name :
Surgical Cautery Equipment

Product Code :
SAMCOMEDICAL1030001A002



Description :

Surgical Cautery Equipment Since the day of our inception, we are engaged in offering Surgical Cautery Equipment (400 W). All the products are designed in compliance with international quality standards, using high-grade raw material and latest technology. Furthermore, we check all the products on various parameters to ensure flawlessness. We also provide customization as per the specifications laid down by the customers. Other information: ½ Feather touch mode selection ½ Delivers precise control of power with various types of current outputs ½ Achieves excellent results in fatty tissue or under water cut ½ Audio visual alarm for separate mode ½ Audio visual patient plate fault alarm ½ MOSFET technology for higher efficiency ½ Independent activation of hand switch or foot switch ½ Monopolar cutting, coagulation, spray coagulation & blended output. & TUR ½ Bipolar coagulation/bipolar cutting ½ Separate generators for monopolar & bipolar

Technical Specification: Monopolar Generator : ½ Frequency: 500 + - 25 Khz ½ Cut Control (Pure): 400 Watts 500 Ohm Load ½ Cut Control (Blend-I): 225 Watts 500 Ohm Load ½ Cut Control (Blend-II): 225 Watts 500 Ohm Load ½ Coagulation Control: 160 Watts 500 Ohm Load ½ Spray Coagulation: 100 Watts 125 Ohm Load Bipolar Generator: ½ Frequency: 500 + - 25 Khz ½ Bipolar Coagulation: 100 Watts 125 Ohm Load ½ Dimension: W 390mm X D 340mm X H 130mm ½ Gross Weight: 27 lbs (approx) ½ Voltage Required: 220V AC 50 Hz ½ Power Consumption: 750 VA Standard Accessories: ½ Double pedal foot switch ½ Patient plate with cord ½ Unipolar handle with cord ½ Hand pencil switch (Disposable) ½ Bipolar forcep with cord ½ Electrode set: 5 No ½ Main Line cord Surgical Cautery Equipment, Surgical Equipment exporters, Surgical Equipment manufacturers

Samco Medical

Website: www.samcomedical.com,

