



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Syringe Needle Destroyer

Product Code :
SAMCOMEDICAL1040001A001

Description :

Syringe Needle Destroyer Our customers can avail from us Syringe Needle Destroyer that is designed in compliance with industry standards & norms. All our products are widely used in hospitals, nursing homes and clinics for destroying needles and syringes immediately after use. Apart from this, we provide these products in different models, sizes and designs at market leading prices to fulfill the demands and needs of customers. Other information: Two slot easy operation for destroying needle & syringe Needle melts electrically & syringe cuts by automatic electrical blades Incinerates the needle instantaneously Cuts the syringe within seconds Low power consumption Can destroy disposable syringe & needle 18G to 28G Collection receptacle to collect the waste Power coating housing Provision of On/Off switch with pilot lamp Shock proof-main circuit breaker is provided Portable, compact & handy Demonstration/Operation : Switch on the device Insert the needle in the needle slot Push it gently downwards till arcing stops Put the syringe in syringe slot Syringe cuts automatically by electrically operated blades Syringe Needle Destroyer, Needle Destroyers exporters, Needle Destroyers manufacturers



Samco Medical

Website: www.samcomedical.com,