



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Three in One Needle & Syringe Destroyer

Product Code :
SAMCOMEDICAL1040001A004



Description :

Three in One Needle & Syringe Destroyer Needle Syringe Destroyer with 3 in 1 operation in single body i.e. burning of needle cutting of syringes and manual cutting of needle are also provided in case of electricity failure. Our machine has strong metal body with powder coating for longer life. Three in One Needle & Syringe Destroyer, Needle Destroyers exporters, Needle Destroyers manufacturers



Samco Medical

Website: www.samcomedical.com,