



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Medical Heating Pad

Product Code :
SAMCOMEDICAL1050001A001

Description :

Medical Heating Pad We have gained expertise in providing our clients with Medical Heating Pad that is available with built-in-snap on-off cut out which helps in maintaining the temperature of the heating pad. Our products are also available with the heating pad to regulate the flow of voltage, which is a dipole switch having option of 0-1-2-3. In addition to this, to keep the pad neat & clean, a washable cover with velcro is provided. Specifications : Size (approx): 9 x 12 inches Wattage (approx): 40 watts Voltage: 220-250 volts (AC) (Operable at low voltage too) Frequency: 50-60 Hz Medical Heating Pad, Heating Pads exporters, Heating Pads manufacturers



Samco Medical

Website: www.samcomedical.com,