



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Baby Incubator

Product Code :
SAMCOMEDICAL130001A004

Description :

Baby Incubator

Technical Specification :

Infant Incubator: Acrylic Single Wall Canopy with six port holes, Front Loading, Four Inlets for administration of O₂, IV lines or Ventilator tubing. Sliding Baby Tray with External jerk free gear based Head Up / Down facility. X-Ray cassette sliding facility. Stainless Steel rust free Heater, blower sub assembly with changeable micro Air filters. Provision of an open shelf in the trolley. Humidity heater and Hygrometer. For Skin / Air / Manual Mode. Control Modes: Skin / Air / Manual Range: Skin control: 34 to 38 C Air control: Room temp. to 39 C Resolution: 0.1 C Accuracy: ± 0.2 C Display: 1in bright skin & air temp.(7 segment LED) Status: On LCD Control Mode & Set temp Timer: Timer / APGAR Timer / Baby Attend Manual mode: Time Programmable Heater output settable Over temp. cut-off Alarm: Auto audio alarm mute, Mute time programmable Safety alarms: ON LCD The mentioned are default limit, can be set or alter as per requirement. Baby Incubator, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,