



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Phototherapy Unit

**Product Code :**  
SAMCOMEDICAL130001A002

**Description :**

Phototherapy Unit

**Technical Specification :**

Phototherapy Unit: Phototherapy Unit with baby bassinet - CFL Phototherapy Stand: 4 Special Blue CFL Tubes. 2 Daylight CFL Tubes. Digital Timer for Lamp Usage Hour and Patient Exposure. Source Cooling Fan. Mirror Coated Reflector. Stand movable on Castors. Infant Care Trolley: Clear acrylic collapsible side panels. Baby tray with transparent slot to incorporate the facility of under surface phototherapy. Mattress with head UP/DOWN facility. X-Ray cassette guide facility. 1 Lower shelves. Movable on castors. Phototherapy Unit, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),